使用说明书
INSTRUCTION MANUAL

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VZUSA

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MEMO

TABLE OF CONTENTS

▪ INTRODUCTION	p. 21
▪ SAFETY PRECAUTIONS	p. 21~22
▪ PRECAUTION BEFORE USE	p. 23
▪ PRECAUTION IN USE	p. 24
▪ PRECAUTION OF MAINTENANCE AND INSPECTION	p. 25
▪ NAME OF PARTS	p. 26~27
▪ EXPLANATION OF FUNCTIONS	p. 28
▪ HOW TO USE	p. 29
▪ TREATMENT PROCEDURE	p. 30
▪ TREATMENT PROCEDURE (Face)	p. 31~32
▪ TREATMENT PROCEDURE (Body)	p. 33~34
▪ PRODUCT SPECIFICATIONS	p. 35
▪ RoHS	p. 36
▪ ASSURANCE PROVISION	p. 37
▪ MEMO	p. 38
▪ WARRANTY CARD	Back cover

INTRODUCTION

- Thank you for purchasing the vzusa SE.

In order to use the device properly and safely, it's necessary to read this instruction manual carefully and follow all safety precautions. Please include this instruction manual if you borrow or give to another person.

- This instruction manual contains the warranty card.

Confirm the information such as " The date of purchase" and " The purchase shop",

SAFETY PRECAUTIONS

Be sure to read this instruction manual carefully before using, and always use properly.

To preserve the risk of injury, loss of life, and damage to any property, the following symbols are used. The marks and meanings of symbols are shown as below.

 DANGER	Failure to observe DANGER signs will result in hazard to personnel is high likely to result in serious injury or death.
 WARNING	Failure to observe WARNING signs could result in hazard to personnel possibly resulting in serious injury or death.
 CAUTION	Failure to observe CAUTION signs could result in injury to personnel and damage to the device and property. ※Property damage includes damage to house, household belongings, domestic animals and pets.

■ MEANINGS OF SYMBOLS

	This symbol indicates prohibition. ※It is indicated a special instruction inside or near a picture.
	This symbol indicates requirement that must be followed.
	This symbol indicates caution.

SAFETY PRECAUTIONS

WARNING

To preserve accidents or discomfort.

Prohibition

■ **Do not use this product with any of the following medical electrical equipment.**

Malfunction of medical electronic equipment causes accidents or discomfort.

- Biological implant such as pacemaker.
- Life support electronic equipment such as artificial heart lung.
- Wearing type medical electronic equipment such as electrocardiograph or drip infusion.

Prohibition

■ **The following persons must not use the product.**

- Persons with heart disease, abnormal blood pressure or malignant tumors.
- Persons with diseases with fever, infectious diseases or acute diseases.
- Persons who are receiving dental treatment.
- Persons with weakened immune systems or who use immune suppressants due to immunity disorders like AIDS etc.
- Persons with hemophilia, bleeding disorders, diabetes or facial neuralgia.
- Persons with paresthesia or dysalgesia.
- Persons with Keloid diathesis, facial melanosis, photosensitivity or pigmentation disorders.
- Persons with capillary dilatation due to long-term steroidal hormone treatments or liver dysfunction.
- Persons who have allergies, atopic dermatitis, skin diseases, sensitive skin or weak skin.
- Persons who have epileptic stroke or discomfort due to flashing lights.
- Persons who have aftereffects of head or ophthalmologic aftereffects due to traffic accidents.
- Persons who have had eye operations (in the case of LASIK, consult with the doctor before use).
- Persons who are pregnant, on period or breast-feeding.
- Children or the person who cannot express their feelings in words.
- Persons under the influence of alcohol or with a hangover.
- Persons who take drugs that cause drowsiness.
- Persons who are extremely tired or do not feel well.
- Persons who have had skin problems such as spots or rashes due to artificial fibers.
- Persons who have or have had skin-inflammation due to cosmetics.
- Persons immediately after having chemical peeling at beauty salons etc.
- Persons who have had skin problems after using beauty equipment.

Prohibition

■ **Do not use on the following parts.**

- Except specific parts.
- Around heats, eyeballs, upper eyelids, lips, and mucous membrane such as inside the mouth.
- Private parts, anus, around the Adam's apple, and ears.
- Parts with morbid spots, moles, and viral warts.
- Tinged and inflamed parts due to excessive sunburn.
- Parts with wounds, eczema, psoriasis, herpes, purulent acnes and/or pimples.
- Parts with pathological spots due to itching, flushing, physical stimulus.
- Parts where surgery was performed. private parts, anus, around the Adam's apple and ears.
- Any parts immediately after hair removal or shaving has been performed.

※ If you are seeing a doctor, consult with the doctor before use.

PRECAUTION BEFORE USE

WARNING

- Do not use this product with other beauty equipment at the same time.
- Do not use this product for any purpose other than beauty and under conditions of malfunction or damage. This product is a treatment equipment specialized for use in beauty. Do not use for other purposes such as therapeutic purposes.
- Do not look at or directly view the LED light.
- Do not use continuously on the same parts.
Problems or symptoms such as skin damage, infection, inflammation, bleeding, pigmentation may be aggravated.
- Do not use serums that contain alcohol ingredients. It may cause skin problems.
- Do not use on pets.
- Do not use near flammable substances such as thinner, benzene, polish removers or sprays.
- Do not use with anything other than the dedicated adapter and USB cord. Do not connect the dedicated adapter and USB cord to other equipment. "The cause of conflagration, electrification, or breakdown."
- Do not damage the USB cord.
- ※ Do not damage, bend more than necessary, pull, twist, place heavy objects on top of, interpose. "Causes of fire, electric shock or damage."
- Do not adopt a star-burst connection. "Causes heat generation or fire."

Prohibition

Prohibition of use
in a wet area

Prohibition of
getting wet

Prohibition of use
with wet hands

Unplug

Compulsion

- Do not use and keep in the bathroom or high humidity places. Do not leave in places where water is easily splashed "on a sink etc."
- Do not immerse in water, pour or spill water (the causes of electric shock, fire or damage).
- Do not touch with wet hands. (the causes of electric shock or injury)
- Make sure to hold the end of the plug when putting in or pulling the adapter. Do not grab the cord to put in or pull the adapter.
- Do not use when the adapter is damaged or the insertion to any outlet is loose.
- Make sure to use only household outlets.
- Insert firmly up to the root of the adapter and the plug.
- Be sure to observe the contents of the instruction manual. Failure to do so may result in some problems.

PRECAUTION IN USE

WARNING

Unplug

- Stop using the product, turn off the power, and pull the adapter from the outlet in the following cases.

Earthquakes or lightning.

Smoke, abnormal heat, odor or distortion of the main body, the adapter or the code.

Prohibition

- During use, this product may slightly heat up but this is not a safety issue.
- Do not bring this product into contact with metals and any other electric appliances such as mobile phone, television.

Compulsion

- Make sure to remove all metal objects, glasses and contact lenses during use.

※Metal objects: Pierced earrings, necklaces, watches etc.

- Make sure to set to level 1 when using for the first time or until you get used to it.
- Do not press the head of product to the skin strongly and never place on the same part of the skin for longer than 3 seconds. This may cause skin troubles.
- During treatments around eyes, use this product with your eyes closed.
- Discontinue use immediately if abnormalities appear on the skin or in the body during use and consult with a doctor.

CAUTION

Precautions of the use of product.

Read the following sections carefully before use.

1

- Perform a special treatment to each part once a week.
- Perform a treatment twice a week in the first and second weeks. From the third week, perform a treatment about once a week on your skin conditions.
- Excessive treatment is not good for the skin.
- Give daily treatments to each part every day or every other day.
- Use this product according to your skin conditions.

2

- As to facial treatments, remove makeup and wash your face to remove impurities such as makeup cosmetics and sebum on your skin and apply the dedicated beauty essence or gel, then use this product.
- Body treatment should be applied after bathing or in other clean skin conditions.

3

- Skin conditions are different depending on your physical conditions.
- Use this product while observing your skin conditions, and discontinue use if abnormalities such as rashes, redness or itchiness appear on skin and consult with dermatologist.

PRECAUTION OF MAINTENANCE AND INSPECTION

WARNING

Prohibition of disassembly

- Do not disassemble, modify, or try to repair the device yourself. For repair, please inquire at dealer.

Unplug

- Unplug the adapter from the outlet when cleaning the main body. (Failure to comply may cause fire or accident.)
- Make sure to turn off the power, and disconnect the AC adapter with the outlet when not using. (Failure to comply may cause fire or damage.)

Prohibition

- Do not adhere dust and pin on AC adapter and plug. (Cause fire, electric shock or damage.)
- Do not drop or hit this product.
- Do not wrap the code around the adapter when storing the device.
- Please keep out of reach of children and do not let them use this product.
- Do not store this product in the following places.
 - Places subject to direct sunlight, high humidity and dust.
 - Places subject to high temperature or temperature below freezing point or inside a car under the blazing sun.
 - Places subject to danger of ignition, rain or wind.
 - Places subject to chemicals, gas generation or excessive vibration.
 - Places being close to heating appliances, appliances generating heat and magnetic sources.

Compulsion

- How to clean

Wipe this product using a soft cloth dampened with water. If wiping with a hard cloth, the product may be damaged or this may cause deterioration in performance.
- ※ Do not use solvents such as alkaline agents, alcohol, thinner or benzene, polish removers, polishing powders, hot water or chemicals.
- ※ Do not leave this product with the beauty essence or gel on it after use.


INSPECTION

When this product is unused for longer than 6 months and is used again, make sure to inspect before use if there are any abnormalities (the lighting state of the LED lights and the presence or absence of cracks).

If used under faulty conditions, this may cause an accident.

NAME OF PARTS

Main Body


※ The images are for reference only. Please refer to the actual product.

Accessories

- AC adapter
- Stand
- Instruction manual with warranty card

NAME OF PARTS


Top cover

- Surface with gold-plated treatment optimized for skin biocompatibility.
- ※ The metal substrate is inherently soft and prone to surface scratches. Please handle with care during use. Minor surface abrasions that may occur will not compromise product integrity and can be safely used as intended.

Equipped with near-infrared (NIR) LED

- This instrument is equipped with near-infrared (NIR) LEDs operating at 830 ± 10 nm wavelength, which activate intermittently after power-on.
- ※ Near-infrared (NIR) light is non-visible with longer wavelengths and enhanced tissue penetration capabilities. Avoid direct eye exposure to the highlighted area, as prolonged exposure may cause ocular damage.

- ※ The white dashed line indicates the location of near-infrared (NIR) LED components.

- ※ Avoid direct eye exposure to the marked area, as prolonged exposure to the light source may cause ocular damage.

EXPLANATION OF FUNCTIONS

Function

Button	Mode	Function
SPA	SPA mode	Through three core mechanisms—thermal energy stimulation, bioelectrical regulation, and physical permeation—it deeply empowers and multi-dimensionally enhances the skin's condition.
UP	ULTRA PULSE mode	It emits high-conductivity pulse waves and employs the "electroporation principle" to enhance the efficacy of substance permeation.
EMS	EMS mode	Output a pulse specialized in working muscles.
MIX	MIX mode	Output insert pulse/EMS/mid-high frequency simultaneously.
VR	VIBRATION mode	Mitigate electric stimulation with the vibration.

LED

	11 Custom-Designed High-Brightness & High-Transparency Red Wave LEDs. For Radiant and Elastic Skin.
	6 Near-Infrared Wave LEDs. Deep-Level Activation to Optimize Skin's Foundation.

- ※ Do not look directly at the LED light.
- ※ Keep your eyes closed during the treatment around eyes.

HOW TO USE

※ When using this device for the first time, remove the protective sticker (transparent sticker) from the electrode.

① Open the cover of the AC adapter socket of the main body and plug the AC adapter into the main unit and the outlet.

② Select the mode

- Press the 〈SPA button〉 to select the SPA mode.
- Press the 〈UP button〉 to select the ULTRA PULSE mode.
- Press the 〈EMS button〉 to select the EMS mode.
- Press the 〈MIX button〉 to select the MIX mode.

※ The mode lamp of the selected mode is lit up in red.

※ Start at level 1.

※ A sound is made once every 10 seconds.

③ Select the level

- This product has 4 types of levels.

※ Select the level by pressing the button of the mode being used.

It is recommended to use at level 1 setting during your first treatment and then adjust to the output level which is comfortable for you.

※ When lowering the level, turn the power of the unit off and then reset the level.

④ VIBRATION selection

- Press the 〈VR button〉 and VIBRATION mode starts.
- VIBRATION mode has two patterns. The order of switching is Short Long OFF.

※ Press the 〈VR button〉 and only VIBRATION mode can be used.

When used alone, the VR mode will automatically stop after 6 minutes of activation.

⑤ Power OFF

- With the activation mode at level 4, press the button again to turn off the power.

※ The SPA mode, EMS mode, and MIX mode will automatically shut down 10 minutes after activation, with a warning tone emitted approximately 30 seconds before shutdown.

※ The UP mode will automatically shut down 6 minutes after activation, with a warning tone emitted approximately 30 seconds before shutdown.

⑥ Disconnect the AC adapter from the main body and outlet.

Close the cover of the AC adapter socket on the main body after use.


TREATMENT PROCEDURE

Before treatment

- 1 Remove accessories such as necklaces or rings and all metal objects that you are wearing.
- 2 Before performing a treatment, remove makeup and wash your face to remove impurities such as makeup cosmetics and sebum on your skin.

Body treatment should be applied after bathing or in other clean skin conditions.

While performing a treatment

- When performing a treatment for full face, do not apply the beauty essence to the face all at once but divide it into four parts to apply as illustrated left and then perform a treatment.

Daily treatment

- EMS mode + MIX mode = 6~10 minutes in total
- Select EMS mode or MIX mode to perform a 6-10 minute treatment.

Special treatment once a week

- SPA mode + UP mode = 10 minutes in total
- 5-minute treatment in SPA mode first, and a 5-minute treatment in ULTRA PULSE mode is effective.


After treatment

- 1 After treatment, gently spread the beauty essence and wipe away the gel left on the skin.
- 2 After treatment, wipe away the product with a soft cloth or tissues.

TREATMENT PROCEDURE (Face)

EMS mode ▪ MIX mode

- Apply the beauty essence on the part of treatment.
- Place this head to the skin in order for approx. 3 seconds as shown in the illustration.
- ①~⑤ Perform each procedure two to three times.

- ※ The treatment procedure in the illustration is reference.
- ※ During the treatment, apply more gel or the beauty essence if the treatment area becomes dry.
When the treatment area becomes dry, it may cause irritation on your skin and the reduction of its effect.
- ※ Please adhere closely 6 elements on skin while treatments.
- ※ During use, wipe the surface of the head using a soft cloth in case the surface of the product's head becomes rough by the adhesion of dirt etc.
- ※ Do not use the product if you are physically unwell.
- ※ Excessive treatment burdens your skin. Do not use multiple times a day on the same area.


WARNING

Do not use around mucosa, upper eyelids, eyeballs, ears, and Adam's apple.

TREATMENT PROCEDURE (Face)

SPA mode ▪ UP mode

- Apply the beauty essence on the part of treatment.
- Move constantly and slowly and do not remain fixed in one place for more than 3 seconds.
- ①~③ Perform each procedure two or three times.

Move along the face line from below the jaw to the rear side of the ear. Then, move from the parotid gland to the shoulder.

Move along the face line from the end of the jaw to the temple. Next, move along the cheek from the nostril to the temple. Finally, move from inner corner of eye, under eye, and to temple.

Move from the bottom to the top while lifting the eyebrows upward to the hairline.

- ※ The treatment procedure in the illustration is reference.
- ※ During the treatment, apply more gel or the beauty essence if the treatment area becomes dry.
When the treatment area becomes dry, it may cause irritation on your skin and the reduction of its effect.
- ※ Please adhere closely 6 elements on skin while treatments.
- ※ During use, in the case that the surface of the product's head becomes rough by the adhesion of dirt etc., wipe the surface of the head using a soft cloth.
- ※ Do not use the product if you are physically unwell.


WARNING

Do not use around mucosa, upper eyelids, eyeballs, ears, and Adam's apple.

TREATMENT PROCEDURE (Body)

Upper arm


- Special treatment: each arm 15 minutes ■ Daily treatment: each arm 15 minutes

- 1 Move in small circles around the root of arm two or three rounds.
- 2 Move straight along the upper arm from above the elbow to the root of arm.

Waist


- Special treatment: 15 minutes ■ Daily treatment: 15 minutes

- 1 Move like drawing small circles around the navel once round.
- 2 Move like drawing small circles from the side of the waist to the inguinal node as if making the curve of the waist.
- 3 Move straight as if flowing from the top to the bottom.

Hip

- Special treatment: 15 minutes ■ Daily treatment: 15 minutes


- 1 Move like drawing small circles from the boundary between the hipline and the thighs to below the waist.
- 2 Move in zigzags around the boundary between the hipline and the thighs.
- 3 Move straight from the boundary between the hipline and the thighs to below the waist.

TREATMENT PROCEDURE (Body)

Leg (front and back)

- Special treatment: each calf 15 minutes
- Daily treatment: each calf 15 minutes


- 1 Move sideways like drawing small circles from above the knee (the back of the knee) to the root of the thighs (the hipline).
- 2 Move linearly like drawing small circles from above the knee (the back of the knee) to the root of the thighs (the hipline).
- 3 Move straight from above the knee (the back of the knee) to the root of the thighs (the hipline).

※ Perform in the same way for the front and the back of the legs.

Calf

- Special treatment: each calf 15minutes
- Daily treatment: each calf 15 minutes

- 1 Move like drawing small circles from the ankle to the back of the knee.
- 2 Move sideways like drawing small circles from the ankle to the back of the knee.
- 3 Move linearly from the ankle to the back of the knee.

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○: Indicates that the content of the hazardous substance in all homogeneous materials of this component is below the limit requirements specified in GB/T 26572-2011.

×: Indicates that the content of the hazardous substance in at least one homogeneous material of this component exceeds the limit requirements specified in GB/T 26572-2011. Currently, there are no alternative technologies or components available, but it complies with relevant national regulations. Meanwhile, our company will continue to strive to reduce the use of hazardous elements through technological improvements.

ASSURANCE PROVISION

- Based on the warranty book records and in accordance with the warranty policy, free repair services will be provided within the warranty period. If a malfunction occurs within the warranty period from the date of purchase, please attach a copy of this warranty book and contact the purchase store or our customer service center for repairs.
- ※ If the serial number, purchase date, purchase store, customer name, address, or phone number is not filled out, this warranty book will be considered invalid. Please ensure all information is completed without omission.
- ※ Accessories (AC adapter, stand and instruction manual with warranty card) are not covered under warranty. Additionally, if the warranty book's serial number sticker is peeled off, the warranty will not apply.
- Even within the warranty period, free repairs will not be provided under the following circumstances:
 - A. The serial number sticker on the warranty book is peeled off.
 - B. The warranty book does not include the serial number, purchase date, purchase store, customer name, address, or phone number.
 - C. Damage or malfunction caused by improper handling by the user, such as dropping or water exposure during transportation or movement.
 - D. Damage or malfunction caused by fire, natural disasters, or abnormal voltage.
 - E. Malfunction caused by connecting to other devices.
 - F. Malfunction caused by using a non-original AC adapter.
 - G. Repairs, adjustments, or modifications performed by non-authorized personnel.
 - H. Failure to present this warranty book.
 - I. Other circumstances.
- This warranty book is only valid for domestic use.
- ※ If purchased overseas, the warranty is only applicable with proof of purchase from an official distributor authorized by our company.
- Regarding Shipping Costs:
 - a. During the warranty period, round-trip shipping costs will be covered by our company. (Valid only within the People's Republic of China.)
 - b. Outside the warranty period, round-trip shipping costs will be borne by the customer.
- This warranty book cannot be reissued if lost. Please keep it safe.

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If there is anything unclear, please feel free to consult at any time.

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